# tRide

Ⓧ	Bedienungsanleitung tRide Rehawagen.....	5
Ⓧ	User manual for tRide stroller .....	15

(1)

(2)

(3)

(4)

(5)

(6)

(7)

(8)

(9)

(10)

(11)

(12)

(13)

(14)


(14)

(15)

(16)

(17)

(18)

(19)

(20)

(21)


(22)

(23)

(24)

(25)

(26)

(27)

(28)

(29)

(30)


(31)

(32)


(33)


(34)



(35)

(36)

(37)

(38)

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### Dear tRide user,

Thank you very much for the trust you have placed in us and for purchasing this product. With your tRide you have acquired an innovative product from the Thomashilfen range.

To ensure that use of the tRide is safe, practical and comfortable for both you and in particular for your child, please read these operating instructions first.

Should you still have questions or problems, please refer to your specialist dealership or contact us directly.

In this context we would like to draw your attention to our website at [www.thomashilfen.de](http://www.thomashilfen.de) where you will also be able to access up-to-date information.

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### Inhalt

	Seite
Intended purpose/indication .....	16
Risks of use and contraindications .....	16
Safety checks and maintenance intervals .....	16
General notes on safety .....	16
Technical specifications .....	17
Setup and assembly .....	17
Parking brake function .....	21
Locking function for the front wheels .....	19
Adjusting the seat tilt angle .....	19
Tip-assit .....	19
Adjustments .....	19
Transporting the tRide .....	20
Product and serial number .....	21
Notes on cleaning, disinfection and care .....	22
Reuse .....	22
Warranty conditions/CE marking .....	22
Maintenance planannual inspection / reuse .....	23

## Intended purpose and places of use

The tRide rehabilitation stroller is suitable for the transfer (at max. walking speed) of children with disabilities by their parents or other approved assistants. It can be used both indoors and outdoors. The seat unit in connection with the subframe may not be used for purposes for which they are not intended (e.g. for transporting other children, shopping, etc.). The tRide rehabilitation stroller is approved for transporting children in motor vehicles/vehicles for people with disabilities in accordance with ISO 7176-19.

## Risks of use and contraindications

Using this product properly will prevent any risks from arising during use. There are no contraindications associated with proper use. The upholstery material used during manufacture of the tRide has been classed as “fire-retardant” in accordance with BS5852 part 1. The extent of resistance to ignition of the materials and assemblies used has been determined in compliance with the standard and has been assessed as uncritical.

## Safety checks and maintenance intervals

- A regular visual check of the operating controls and securing screws should be carried out by a trained specialist every month.
- Please remember to re-secure the screws, star handles and locking levers (hip angle adjustment) following each adjustment.

## Please observe the following notes on safety

- Before using the tRide for the first time, read the operating instructions carefully – or have them read aloud to you if you have difficulty reading. If you have lost the operating instructions you can download them at any time from the download page of our website at [www.thomashilfen.de](http://www.thomashilfen.de)
- We recommend that you readjust the settings of the tRide at least every three months to suit the actual size of your child's body. This should be carried out with assistance from your therapist/orthopedic technician where appropriate.
- **CAUTION:** The gap between the chassis and the brake actuating rod presents a potential risk of entrapment of clothes and limbs in the foot area.
- **CAUTION:** The gap between the seat and the chassis frame presents a potential risk of entrapment of clothes and limbs in the seat area.
- **CAUTION:** The effect of solar radiation, electric heaters or other heat sources can cause the seat/rehabilitation stroller to overheat. This presents a potential risk of burning to the back and legs upon contact! Always ensure therefore, that the seat/rehabilitation stroller is not subjected to such sources of heat during use or storage and only use it after it has cooled sufficiently.
- **CAUTION:** Wet or flat tires can affect the performance of the brakes. During a stop or longer pause you should secure the rehabilitation stroller against unintended rolling by applying the foot brake.
- Only lift the rehabilitation stroller by gripping securely welded or bolted components. (Front frame tube above the front wheels, rear axle, push handles/push bar).
- Ensure that you do not injure yourself when folding and unfolding the rehabilitation stroller.
- Please do not overload your tRide rehabilitation stroller and observe its maximum load limitation (see technical specifications).
- We recommend that you always secure the child using the 2-point pelvic belt or with another type of restraining strap (available as an accessory) and that you never leave the child unsupervised in the rehabilitation stroller.
- Never allow your child to get into or out of the rehabilitation stroller unsupervised. If the footrest is overloaded there is a risk of slipping or tipping over.
- High-visibility reflective clothing will make you and your child easier to see for other road users at night.
- Heavy bags or shopping nets attached to the rehabilitation stroller will increase the risk of it tipping over. Instead use the basket or tray (available as accessories) on the subframe.
- When using the tRide please pay particular attention to the recommendations in the operating manual regarding folding and unfolding the stroller.
- Please keep your children away from the packaging – there is a risk of suffocation.
- Only one child may be placed in the rehabilitation stroller.
- Market experience to-date does not indicate any dangerous combinations of parts/accessories of the tRide rehabilitation stroller.

## Technical specifications

	tRide 1	tRide 2
Seat depth	17 - 28 cm / 6.7 - 11"	24 - 35 cm / 9.4 - 13.8"
Seat width	20 - 30 cm / 7.9 - 11.8"	26 - 35 cm / 10.2 - 13.8"
Back height*	26 - 42 cm / 10.2 - 16.5"	32 - 52 cm / 12.6 - 20.5"
Lower leg length	17 - 30 cm / 6.7 - 11.8"	25 - 36 cm / 9.8 - 14.2"
Hip angle	90° bis 140°	90° bis 140°
Knee angle	90° bis 180°	90° bis 180°
Foot angle	78° bis 100°	78° bis 100°
Height adjustment arm support	10 - 14 cm / 3.9 - 5.5"	12 - 18 cm / 4.7 - 7.1"
Seat tilt - in-space	-10° bis +35°	-10° bis +35°
Size of arm support(w x h)	21 x 10 cm / 8.3 x 3.9"	29 x 12 cm / 11.4 x 4.7"
Seat height	46 cm / 18.1"	47 cm / 18.5"
Size of footrest (w x d)	29 x 18 cm / 11.4 x 7.1"	33 x 21 cm / 13 x 8.3"
Overall dimensions of seat unit (l x w x h)	85 x 41 x 57 cm 33.5 x 16.1 x 22.4"	98 x 46 x 46 cm 38.6 x 18.1 x 18.1"
Folded dimensions of seat unit (l x w x h)	63 x 36 x 38 cm 24.8 x 14.2 x 15"	76 x 46 x 50 cm 29.2 x 18.1 x 19.7"
Weight (seat)	8,8 kg / 19.4 lbs	10,8 kg / 23.8 lbs
Push bar height	76 - 117 cm / 29.9 - 46.1"	81 - 120 cm / 31.9 - 47.2"
Overall dimensions of chassis in running order (l x w x h)	75 x 60 x 77 cm / 29.5 x 23.6 x 30.3"	81 x 73 x 79 cm 31.9 x 28.7 x 31.1"
Folded dimensions of chassis (l x w x h)	79 x 60 x 38 cm 31.1 x 23.6 x 15"	83 x 67 x 38 cm 32.7 x 26.4 x 15"
Overall dimensions of the tRide rehabilitation stroller, ready for use (l x w x h)	90 x 60 x 101 cm 35.4 x 23.6 x 39.8"	110 x 67 x 99 cm 43.3 x 26.4 x 39"
Folded dimensions of the tRide rehabilitation stroller (l x w x h)	81 x 60 x 63 cm 31.9 x 23.6 x 24.8"	79 x 64 x 72 cm 31.1 x 25.2 x 28.4"
Wheel size (front (swivel wheels))	7,5"	7,5"
Wheel size (rear)	10"	10"
Weight (chassis without seat unit)	8 kg / 17.6 lbs	9 kg / 19.8 lbs
Carrying capacity (seat)	35 kg / 77.2 lbs	40 kg / 88.2 lbs
Carrying capacity (chassis)	45 kg / 99.2 lbs	50 kg / 110.2 lbs

\*without headrest

## Setup and assembly

Settings/adjustments may only be undertaken by specialists that have received the requisite training or instruction (e.g. a medical product consultant at a medical appliance specialist dealership).

Fitting/removal of accessories is carried out using commercially available standard tools. Separate assembly and operating instructions are supplied with each accessory. If you do not have these instructions you can download them at any time from the download page of our website at [www.thomashilfen.de](http://www.thomashilfen.de)

### Unfolding

Hold the rehabilitation stroller by its push handle, place your foot on the rear axle and then pull the handle upwards (until the push handle audibly clicks into place) (fig. 1). If necessary, tilt the seat to the horizontal position (see “Adjusting the seat tilt angle”). Then bring the backrest into the upright position (release the locking lever on the back of the seat unit). The footrest must be secured at the toggle joints. To do this, place the ratchet joints one above the other and secure with the screws (from the inside), the washers (from the outside) and the wing nuts (from the outside).

### Folding

Bring the backrest into the prone position, release the toggle joints and swing the footrest towards the seat unit. Tilt the seat to the horizontal position (see “Adjusting the seat tilt angle”) and then release the locking levers (pull upwards) located on both sides of the subframe while at the same time pushing the push handle downwards and inwards toward the frame (fig. 2).

**CAUTION:** Hold the frame securely – do not allow it to fall. Do not apply the foot brake; otherwise there is a risk of deforming the wheels.

**Note on folding:** Ensure that the seat is always in the horizontal position when folding your tRide.

### Fitting the seat unit

The lever on the headrest frame serves on the one hand to separate the seat unit from the subframe and on the other to adjust the seat tilt angle. Hold down the lever on the headrest frame of the seat unit (fig. 3).

Insert the seat unit into the mounting adapter on the chassis and ensure that the pins on the side of the seat unit engage in the recesses designed for this purpose. Bring the seat into the upright position until you hear it securely click into place. Then check that the seat is properly locked in-place by pulling forwards on it (fig. 4).

### Removing the seat unit

The seat can be removed from the subframe without the use of tools. Hold down the lever on the headrest and tilt the seat unit backwards. The seat unit can then be removed by lifting it diagonally backwards and upwards out of the mounting adapter on the chassis. It is recommended that you first lock the foot brake in the on position.

### Fitting the footrest

The footrest must be secured at the toggle joints. To do this, place the ratchet joints one above the other and secure with the screws (lever the inside), the washers (lever the outside) and the wing nuts (lever the outside) (fig. 5).

### Lateral backrest upholstery

Place the lateral backrest upholstery over the lateral tube and use Velcro to fasten the upholstery to the curved lateral parts of the back section (fig. 6). Completely close the Velcro – including over the knurled screw. Also use Velcro to secure it to the back section (fig. 7).

### Locators for weather canopy

The tRide rehabilitation stroller is fitted with locators to fit a weather canopy as standard (fig. 8)

### Seat upholstery

Use Velcro to secure the seat upholstery to the seat. The two lower Velcro straps should be secured around the footrest tubes (fig. 9).

### Footrest tube padding

Use Velcro to secure the footrest tube padding around both the left and right tubes beneath the toggle joints.

### Mounting the belts

Pass the two upper belts depending on the size of the child to either both inner or both outer loops of the back upholstery and attach them to the rear of the backrest with the cam buckles. (Fig. 10)

The two bottom belt ends can be attached with the cam buckles to either the back plate or to the side of the plastic element. (Fig. 11-13) Guide the belts left and right next to the back upholstery to the corresponding cam buckle and fix them.

**Tip:** The belt length can also be adjusted from the front. Pull the D-ring to shorten the belt.

### Basket

Fasten the basket (with its opening to the rear) to the four crescent rings on the frame using the spring clips (fig. 14). The basket may remain attached when folding the tRide. (Maximum loading of the basket: 5 kg).

## Parking brake function

Please remove the cable tie that locks the foot brake.

**Lock:** Place your foot on the brake bar and press downwards.

**Release:** Place your foot underneath the brake bar and slide it upwards.

**CAUTION:** The brake release force is higher than usual in order to provide improved braking performance.

## Locking function for the front wheels

**Lock:** Straighten the front wheels and lock them in-place using the locking pin (audibly clicks into place).

**Unlock:** Pull out the locking pin and fix it using a half-turn (fig. 15)

## Adjusting the seat tilt angle

**CAUTION:** When adjusting the seat tilt angle always hold the back of the rehabilitation stroller to prevent the child from being jolted when changing position. We recommend that you make this adjustment with the child not sitting in the seat unit. If the adjustment is made when the child is sitting in the seat the child must always be held or supported by a second person.

The seat tilt angle is adjusted using the lever on the headrest frame. Pulling and holding the release lever releases the seat tilt mechanism which can then be set to a range of positions. These positions are evident from the hole pattern on the seat mount. Once you have positioned the seat as desired, release the lever and the seat tilt angle will then automatically lock into place. Check once again that the seat has locked into place properly by pressing on the seat's backrest. The seat tilt mechanism must be locked into the corresponding locking hole on both sides of the tRide (fig. 16)

**CAUTION:** During day-to-day use always ensure that the release lever remains unobstructed and that it is not accidentally released or locked (by pressing against it or through it coming into contact with other objects). There is otherwise an increased risk of tipping over/injury to the child.

Active seating position: The seat tilt adjustment of the tRide can be set to an active seating position of  $-10^\circ$ . A non-locking spring-loaded pin prevents this position from being accidentally selected. To select the active seating position, pull out the spring-loaded pin on the seat mounting while adjusting the seat tilt angle (fig. 33)

Once the user leaves the active seating position the spring-loaded pin will once again block the seat tilt mechanism when the  $0^\circ$  position is exceeded.

### Tip-assist (accessory)

Place your foot on the tip-assist then – with both hands on the rehabilitation stroller – tip as required (fig. 17).

## Range of adjustments

### Seat width

Loosen the four star grips (fig. 18) beneath the seat and slide out the side panels as necessary (fig. 19). If the seat depth is set to its minimum, the maximum setting for the seat width is governed by the available space between the lateral backrest upholstery pads. The seat upholstery does not need to be removed to adjust the seat width. Once you have made the adjustment, re-tighten the star grips.

### Seat depth

Loosen all four star grips (that are also used to adjust the seat width) beneath the seat and pull out the seat section towards you (fig. 20). The seat upholstery does not need to be removed to adjust the seat depth. Once you have made the adjustment re-tighten the star grips.

**Tip:** For easy adjustment of the seat depth grip the front edge of the seat with one hand and press against the bottom part of the backrest (center) with the other.

### Arm support

The arm supports are adjustable both in height and depth. To adjust the arm supports loosen both wing nuts, set the desired height and depth and retighten the nuts (fig. 21).

### Back height

Slightly release the lateral backrest upholstery from the seat. Loosen the two knurled screws at the top of the lateral tubes (fig. 22). Now the back section can be infinitely adjusted in height until the point where it locks into place (the spring-loaded pin snaps audibly into place). Use Velcro to re-attach the lateral backrest upholstery.

To reduce the back height, press the spring-loaded pin beneath the knurled screw on the right side and slide the back section downwards until it reaches the desired height (fig. 23).

## Headrest

Insert the headrest tube into the locator on top of the headrest bracket and tighten the wing screw (fig. 24). The headrest can be infinitely adjusted in height using the knurled screws on the rear of the backrest (fig. 25). The headrest angle can be changed by loosening the wing screws on the headrest bracket (fig. 26).

Furthermore, the depth of the headrest can be adjusted by loosening the wing screw on top of the headrest bracket (fig. 24). The angle of the headrest surface can also be adjusted both vertically and horizontally by loosening the screws of the ball joint (fig. 27). Ensure that all screws are retightened after the adjustments have been completed.

## Back angle adjustment

Release the locking lever and pull out the safety locking pin located beneath the lever (fig. 28). The backrest can now be set to the desired angle. It must be ensured that the locking pin snaps back into place. Finally, return the locking lever to the locked position once the adjustment has been made.

**Note:** The seat units has got a number of hole on the backrest into which the safety locking pin snaps.

When transporting the rehabilitation stroller in a vehicle, the safety locking pin must always be located in the lowest hole (90° seat angle).

## Lower leg length

Loosen the threaded pins at the side of the foot rest using the Allen wrench (fig. 29). The foot rest can now be can be infinitely adjusted. Finally, retighten the threaded pins once the adjustment has been made.

## Foot angle

Loosen the four screws on the right and left on the underside of the foot rest. Adjust to the desired foot angle and retighten the screws (fig. 30).

## Knee angle

Set the knee angle by loosening the wing nuts located on both toggle joints (fig. 31).

**Tip:** The foot rest can be flipped up to facilitate transfer of the child.

## Transporting the tRide

The tRide can be folded for transport in the luggage area of a vehicle.

The tRide has been successfully impact tested against international standard ISO 7176/19 and ANSI/RESNA WC/Vol.1 – Section 19 6/22/99 Draft and may be used as a vehicle seat.

We recommend that the transport of children with a body mass of less than 22kg transfer to a Child Restraint System (sometimes referred to as an infant safety seat, a child restraint system, a restraint car seat, or ambiguously as car seats) in compliance with ECE Regulation 44/04; In order to provide the best possible protection for the child.

If transfer to a Child Restraint System is not possible, due to the needs of the occupant e.g. a need for continual postural support, occupant discomfort, anxiety or a dependency on the medical device, then the decision to carry the occupant in their wheelchair or buggy should involve an individual risk assessment process. This risk assessment process should involve a therapist who would be able to support reasons not to transfer the occupant to a Child Restraint System.

**IMPORTANT:** When used as a vehicle seat the following instructions must be followed.

The tRide is designed to be forward-facing when used as a seat in a motor vehicle. The occupant should transfer to the vehicle seat and use the vehicle-installed restraint system whenever it is feasible!

In order to fasten the chassis at the floor of the vehicle, only use belt systems which meet the standard ISO 10542.

The 4 tie-down securement points are indicated with the following symbol:

**CAUTION:** When the stroller is being used as a seat in a motor vehicle the locking pin at the tilt-in-space mechanism must engage in a position that inclines the seat backwards. (fig. 32)

## Preparing the tRide for transport:

1. Remove the following (if fitted); pommel, play tray, rain hood, sun canopy, grip rail, accessory bag and shopping basket contents.
2. Position the stroller correctly in the vehicle over the floor track with the front of the stroller facing towards the front of the vehicle (in the direction of travel). Put the brakes on.
3. Make sure the hip angle adjustment is locked. The locking pin at the backrest must be engaged (see marking).

## Securing the tRide in a Vehicle:

1. The stroller must be secured in a vehicle using a 4-point tie down system that complies with either ISO 10542 Part 2, SAE 2249 or DIN 75078/2. (Such restraints are manufactured by Unwin Safety Systems and Q'Straint)
2. Vehicle operators and/or assistants should be fully trained in the use of wheelchair tie-down and occupant restraint systems.
3. The two front tie-downs must be attached left and right side to the steel ring located on each front castor mounting boss (fig. 34).
4. The two rear tie-downs must be attached to the two securement points located on each side of the pushchair frame (fig. 35).

## Fitting the occupant restraint system:

1. The stroller passenger must be restrained using an approved vehicle safety belt. Thomashilfen recommends the use of either a split reel double inertia seat belt or a 3-point occupant restraint system. (Such restraints are manufactured by Unwin Safety Systems and Q'Straint)
2. Vehicle operators and/or assistants should be fully trained in the use of wheelchair tie-downs and occupant restraint systems.
3. The occupant lap belt must be routed low over the pelvis, avoiding contact with the abdomen of the passenger. Lap belt routing should be inside of the occupant side pad on both left and right sides of the pushchair (fig. 36).
4. The lap belt should not be held away from body by rehabilitation stroller components or parts, such as side pads or wheels (fig. 37).

**NOTE:** The belt webbing should not be twisted when in use.

5. The angle of the pelvic belt should be within the optimal zone of 30° to 75° to the horizontal (fig. 38).
6. The occupant lap belt should be tightened to make a snug fit without causing discomfort to the passenger.
7. If the pushchair is fitted with a harness, this should only be fastened after the occupant restraint system.
8. The shoulder belt of the occupant restraint system must be fitted to the side wall of the vehicle at a height that ensures the seat belt webbing lies in the middle of the occupants shoulder.

**NOTE:** To meet current regulations only a shoulder belt fitted to the side wall of the vehicle should be used to correctly restrain the occupant.

**Product/serial number** (references the year of manufacture)

To ensure you receive the correct replacement parts or for potential technical product recalls we will always require you to state the serial/batch number together with the exact item number. The serial/batch number will allow you to find out the year of manufacture of your tRide from us. The serial/batch number and the item number of the tRide can be found on the barcode label located beneath the seat.

**Notes on cleaning, disinfection and care**

**Subframe**

The tRide subframe is manufactured from high quality steel and aluminum tube. It is finished with a durable powder coating that can be cleaned using household detergent (neutral soap).

**Upholstery**

All of the upholstery can easily be removed. The 100% PES material and the yarn interlacing are made of high-quality fabric are washable at up to 40°C. The materials are printed with non-toxic dyes that are free of heavy metals. We draw your attention to the fact that even high-quality materials may fade as a result of sustained and intensive exposure to solar radiation or due to frequent washing. The upholstery is fire retardant (BS 5852 part 1).

**Velcro fasteners**

In order to maintain the functionality of the Velcro fasteners they should be brushed down from time to time. Always close the Velcro fasteners before washing.

**Disinfection**

Clean the upholstery and grips before disinfecting. Wipe down all parts of the product with a cloth dampened with disinfectant. Only water-based disinfecting agents should be used. The manufacturer’s instructions for use should be observed when disinfecting.

**Reuse**

The tRide is suitable for reuse. However, this rehabilitation product is subject to exceptional stresses. We have calculated the product’s service life – when used properly and in compliance with the service and maintenance schedule – to be 5 years, on the basis of market observations and the state of the art. Periods during which the product is in storage at specialist dealerships are excluded.

With corresponding care and maintenance this product will also be reliable in operation for significantly longer than this period.

Before reuse the product must always be carefully cleaned and disinfected. The condition of the product must be checked by trained specialist personnel for wear and damage. Any problems with the product must be rectified in order to ensure that the rehabilitation stroller can continue to be used safely.

If the condition of the rehabilitation stroller is unsuitable for reuse it can be recycled at your local recycling facility.

**Guarantee conditions / CE mark**

Thomashilfen offers you a guarantee of 3 years for all frame parts, commencing on the date of purchase. This guarantee incorporates all claims which influence functionality. Damages due to improper use (e.g. overloading), and natural wear and tear are excluded.

Swiftly meets the requirements of the European standards EN 12182 and 12183 and is provided with a CE marking. Enjoy using your Swiftly rehabilitation buggy!

**Maintenance plan-annual inspection / reuse**

A user manual should always be delivered with the product. If the user manual is missing please contact the manufacturer to request a new user manual. Please make yourself familiar with the product. If you don’t know the product and its functions please read the user manual carefully before inspection. Please clean the product before inspection. Pay attention to all caring and all product specific instructions mentioned in the user manual.

Pos.	Part <input checked="" type="checkbox"/> = checked	(Service / Reinstatement)					
		Function	Damage	Distortion	Replace		
1	Seat	S/R	S/R	S/R			
2	Angle adjustment	S/R	S/R	S/R			
3	Seat upholstery		S				R
4	Missing caps and plugs						S/R
5	Back	S/R	S/R	S/R			
6	Angle adjustment	S/R	S/R	S/R			
7	Height adjustment	S/R	S/R	S/R			
8	Back upholstery		S				R
9	Headrest	S/R	S/R	S/R			
10	Headrest upholstery		S				R
11	Push bar	S/R	S/R	S/R			
12	Push bar coating		S				R
13	Chassis	S/R	S/R	S/R			
14	Tyres	S	S	S			R
15	Parking brake	S/R	S/R	S/R			
Options							
16	Lateral trunk supports	S/R	S/R	S/R			
17	Lateral trunk support upholstery		S				R
18	Arm supports	S/R	S/R	S/R			
19	Arm supports upholstery		S				R
20	Abduction block	S/R	S/R	S/R			
21	Abduction block upholstery		S				R
22	Tray	S/R	S/R	S/R			
23	Harness	S/R	S/R	S/R			
24	Other accessories	S/R	S/R	S/R			
Service done by: :				Date:			

Servicing and repairs to the rehabilitation stroller may only be carried out by a specialist dealership. The user of this rehabilitation aid must ensure that the prescribed inspections/maintenance are/is performed regularly and at the correct intervals. If damage is evident the user is under a duty to actively inform and to have the damage repaired immediately by specialist personnel. Only original spare parts may be used for repairs (the list of replacement/spare parts can be found on the download page of our website at [www.thomashilfen.de](http://www.thomashilfen.de)).

The removal and fitting of spare parts may only be carried out by specialist personnel! If it becomes necessary to return the tRide to your specialist dealership, please ensure that it is properly and securely packed for transportation purposes.

Member of the International  
Support Association for the  
Rehabilitation of Children and Youths e.V.

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Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG