

Defender ReHA CAD Order Form 01.2021

Date	Buyer	PO#	Mark for
Address			
Contact	Phone	Email	

Step 1 - Choose Model

Model 888 - CAR Seat/Back with height-adjustable headrest, LATCH equipped, 5-point safety/ positioning belt with shoulder pads, chest clip, small crotch pad, set of removable head, trunk and hip supports, cup holder, two-tone colour - black-charcoal.

Model 889 - BOOSTER Seat/Back with height adjustable headrest, 5-point positioning belt with shoulder pads, chest clip, large crotch pad, set of removable head, trunk and hip supports, cup holder, two-tone colour- black-charcoal.

888	Car Seat	\$1499
-----	----------	--------

887	Booster Seat	\$1499
820-1	Turning Plate	\$899

Step 2 - #888 Car Seat Options

FOOTREST OPTIONS :		Code	Price
Footrest Adapter + Base plate	Ea	821	\$550
Footrest, with short footplate (depth 4.7") (tilt/ angle adjust.)	Ea	806	\$410
Footrest, with long footplate (depth 7.1") (tilt/ angle adjust.)	Ea	805	\$410
Footrest pad (for 806), height 4", black	Ea	834	\$125
POSITIONING SUPPORTS:			
Abduction block, black	Ea	833	\$160
Extended Crotch Strap *	Ea	814	N/C
Seat Wedge, below, black, 10 degree	Ea	810	\$250
Table, black	Ea	831	\$415

Step 2 - #889 Booster Seat Options

FOOTREST OPTIONS :		Code	Price
Footrest Adapter + Base plate	Ea	821	\$550
Footrest Short (depth 4.7") (tilt/ angle adjust.)	Ea	806	\$410
Footrest Long (depth 7.1") (tilt/ angle adjust.)	Ea	805	\$410
Footrest Pad (for 806), height 4", black	Ea	834	\$125
POSITIONING SUPPORTS:			
Abduction Block, black	Ea	833	\$160
Extended Crotch Strap	Ea	814	N/C
Seat Wedge, below, black, 10 degree	Ea	810	\$250
Table, black	Ea	831	\$415
OTHER OPTIONS:			
Booster Latch Connection Kit for turning plate	Ea	818	\$190
Front Stand	Ea	817	\$250
Stabilizing Bow	Ea	819	\$125

* Extended Crotch strap is only to be used when the car seat is converted to Booster mode, it is not safe to use in Car Seat Mode.

Technical Specifications	
Child weight	888 - Car seat - 22 - 65lbs 887 - Booster - 40 - 110lbs
Child height	888 - Car seat 27 - 57" 887 - Booster 34 - 57"
Seat weight	22 lbs
Seat depth	13.5"
Seat width	8.5 - 12"
Seat to head	23.5 - 30.50"
Shoulder height	10-19.5"
Lower leg length	7" with footrest pad (#834) 11 - 17" Standard (#805/806)



The Defender is a Universal Restraint System for disabled children.

It is approved in accordance with the Canada Motor Vehicle Restraint System and Booster Seat Safety Regulations, CMVSS213, for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System for this age group.

Only suitable for use in vehicles with 3-point (lap & shoulder) seat belts. Suitable only for forward facing installations.

The Car seat is for children weighing between 10 - 29.5kg (22 - 65 lbs.), whose height is between 68.5 - 144.7cm (27 - 57") and who are over 12 months old.

The Booster seat is for children weighing between 18.14 - 49.8kg (40 - 110 lbs.) and whose height is between between 86.3 - 144.7cm (34 - 57").

