

Swiftly CAD Order Form 10 2019

Date:	Contact:	PO#	Marked For:
Customer Name:			
Address:			
Confirmation Email To (Buyer):		Phone:	Fax:


Swiftly
 Swiftly
Seat Size 1 - \$2500
Seat Size 2 - \$3195
Description

680621 - Seat	680821	Folding frame with integrated/contoured seat, fixed tilt, adj. seat depth from underneath, footrest height adjustable, back recline, suitable for early intervention use.
680622 - Back	680822	
6830 - Bus Tie downs	6830	
66295 - H-belt	66295	

	Lower Leg	Seat Depth	Seat Width	Back Height	Hip Angle
Size 1	6.25 - 13"	8.25 - 11.25"	13.5"	24.4"	90 - 115 Degrees
Size 2	8.3 - 14.6"	14.8"	15.6"	27.6"	90 - 100 Degrees

Seat Size Chart

Lower Body Positioning

	Unit	Small	Medium	Large	Price
Foot Straps with Cam Buckle	Pair		6648		\$45
Ankle Huggers, Bodypoint	Pair		213151		\$250
Footrest Pad, Height 1 3/4	Each			10023	\$115

Mid Body Positioning

Abduction Block, Padded	Each	6620		6720	\$185
Hip Supports 2 Point, Padded SR Buckle	Each	8420		8470	\$90
Hip Supports Soft Pelvic Harness	Each	6614		6714	\$195
Crotch Strap (fits 6629, 8420, 8470, 8421, 8471)	Each	6341		6341	\$30
Tray, Height Adj, Clear with Picture Window	Each	68221		67281	\$545
Mounting Hardware Tray *Req'd For 6628/6728	Each	68222		67282	included

Upper Body Positioning

Chest Butterfly Style Neoprene	Each	6385		6585	\$215
Chest Vest Style	Each	6631	6632		\$185
H-style, w/integrated hip belt			66295		\$130
H-style, posture support, AEL	Each		14813		\$225
Pivot shoulder harness, dynamic,	Each		214270		\$225
Pivot shoulder harness	Each		214250		\$200
Trunk Support, Rigid Lateral	Pair	6634		6734	\$300
Trunk Support, Flex Lateral w/ Chest Strap	Pair		6633	6633	\$420
Contoured pillow	Each	6824		6837	\$215

General Accessories

Seat minimizer	Each	6834		6836	\$400
Pneumatic wheels, full set	Each	6646		6671	\$275
Cargo basket	Each	6813		6835	\$95
Frame Padding	Pair			6622	\$155
Grab Bar, Padded Black	Each	6833		6721	\$140
Tip Assist	Each			6829	\$95
Straps, Secure folded seat/frame	Each			667011	\$35
Straps, Wrist "Runaway" Tether	Each			667111	\$35

Weather Protection

Canopy, Sun/Rain with Window Black	Each			663761	\$375
Lower Body Wind Breaker, black for 663761	Each			661852	\$130
Canopy, Rain Only, Clear	Each			6638	\$250
Rain Shield, Clear for Canopy 663761/637811	Each			6739	\$150
Rain cape, Grey	Each	6340		6540	\$200
Sleeping bag, Summer Weight, Black	Each	6619		6719	\$385
Sleeping Bag, Winter Weight, Black	Each	6618		6718	\$400
Umbrella, Angle Adjustable	Each			6594	\$115

Swiftly	Size 1	Size 2
Frame capacity	77lbs	110lbs
Unit weight	27lbs	40lbs
Max seat tilt	15 degrees	20 degrees
Back recline	90-115 degrees	90-100 degrees
Upholstery	Trevira-CS	Trevira-CS
Floor - handle	28-45"	32.3 - 48.8"
Floor to seat	19"	
Max width	24"	26.4
Fold length	28.75"	31"
Fold height	15"	19.3"
Brakes	Foot operated	Foot operated
Seat direction	Front facing	Front facing
Front Wheels	7.5" solid	7.1" solid
Rear Wheels	10" solid	9.8" solid

Seat Minimizer Specs

	Size 1	Size 2
Seat Depth	6.9-11.2"	11-14.8"
Seat Width	9.1"	11.6"
Footrest Pad Height	1.8"	2"

