

ThevoTwist CAD Order Form 1.0

Date:	Contact:	PO#	Marked For:
Customer Name:			
Address:			
Confirmation Email To (Buyer):		Phone:	Fax:

Step 1 - Seat Sizes (required options indicated by * symbol)



	ThevoTwist Seat Size 1 \$2,500		ThevoTwist Seat Size 2 \$2,600		Description
	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	
Seat*	668821	668811	678821	678811	Individual body adapted, integrated seat/back with MIS-Microstimulation®; 145 degree recline gas spring operated; growth adjustable seat upholstery; articulated footrest height and angle adjustable swing away head rest bracket; front/rear facing seat; ventilated/washable upholstery in 2 colors.
Back*	668822	668812	678822	678812	
Recline*	668823	668813	678823	678813	
Hip Sup.	668824	668814	678824	678814	
Calf Strap*	6635	6635	6735	6735	



6387/6587 Standard* 6598 Neck Support* 6599 Shell*

*Headrest is a required option and included in seat base price

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Back Height #	Max Height	Hip Angle
Size 1	7 - 11.5"	6.75 - 12"	7.75 - 11"	18.25 - 26.75"	49"	90 - 133 Degrees
Size 2	8.5 - 14.5"	9.25 - 15"	10 - 13.5"	25.25 - 30.75"	60"	90 - 145 Degrees

Step 2 - Frame Options (Required options are indicated by * symbol)

A-chassis - Folding aluminum frame, tilt-in-space, dual spring suspension, basket, tilt-in-space lever at push bar

	Size 1	Size 2
A-chassis	6202 <input type="checkbox"/> \$2,595	6402 <input type="checkbox"/> \$2,695
Castor Pin Lock*	67830 <input type="checkbox"/> included	67830 <input type="checkbox"/> included
Bus Tie-Down Loop Kit*	67839 <input type="checkbox"/> included	67839 <input type="checkbox"/> included
Pneumatic Wheels, Full Set	6646 <input type="checkbox"/> \$260	6646 <input type="checkbox"/> \$260
Frame Padding, Reversible	6622 <input type="checkbox"/> \$155	6622 <input type="checkbox"/> \$155
Umbrella with A-chassis	6594 <input type="checkbox"/> \$115	6594 <input type="checkbox"/> \$115
Tip Assist	6750 <input type="checkbox"/> \$95	6750 <input type="checkbox"/> \$95
Cargo Basket	included	included



ThevoTwist Technical Specs **

	Size 1	Size 2
Seat Capacity	77 lbs	88 lbs
Frame Capacity	99 lbs	110 lbs
Unit Weight	36 lbs	42 lbs
Max Seat tilt	31"	
Upholstery	Trevira-CS	
Foot Angle	78" - 100"	
Size of Foot Rest	11.5" x 7"	13" x 8"
Floor to Seat	19"	19.5"
Max Width	23"	27"
Fold Length	32"	35.5"
Fold Height	18.5"	18.5"
Brakes	Foot Operated	
Seat Direction	Front and Rear Facing	
Front Wheels	7.5" solid	
Rear Wheels	10" solid	

Q-chassis - Hi-low indoor base/feeder/therapy frame, hght/tilt adj. w/ pedal gas spring, lock rear wheels

	Size 1	Size 2
Q-chassis	6666 <input type="checkbox"/> \$2,365	6666 <input type="checkbox"/> \$2,365
Push Bar for Q-chassis	6660 <input type="checkbox"/> \$370	6660 <input type="checkbox"/> \$370

Q-chassis Technical Specs**

Capacity	88 LBS
Seat Tilt	-7 to +31°
Floor to Seat	11 to 28"
Seat Direction	Front Facing



Z-chassis - Indoor base/feeder/therapy frame, manual height/tilt adjustable lockable rear wheels

	Size 1	Size 2
Z-chassis	6691 \$1,500	6691 \$1,500

Z-chassis Technical Specs***

Capacity	88 LBS
Seat Tilt	-15 to +35°
Floor to Seat	15 - 23"
Seat Direction	Front Facing



ThevoTwist CAD Order Form 1.0

Step 3 - Positioning and General Options

Lower Body Positioning	Unit	Small	Medium	Large	Price	Order
Foot Straps with Cam Buckle	Pair		6648 <input type="checkbox"/>		\$45	
Ankle Huggers, Bodypoint®	Pair		213151 <input type="checkbox"/>		\$250	
Footrest Pad, Height 1 3/4	Each	10622 <input type="checkbox"/>		10023 <input type="checkbox"/>	\$115	

Mid Body Positioning

Abduction Block, Padded	Each	6620 <input type="checkbox"/>		6720 <input type="checkbox"/>	\$185	
Hip Supports 2 Point, Padded SR Buckle	Each	8420 <input type="checkbox"/>		8470 <input type="checkbox"/>	\$90	
Hip Supports 4 point, Padded SR Buckle	Each	8421 <input type="checkbox"/>		8471 <input type="checkbox"/>	\$245	
Hip Supports Soft Pelvic Harness	Each	6614 <input type="checkbox"/>		6714 <input type="checkbox"/>	\$195	
Crotch Strap (fits 8420, 8470, 8421, 8471)	Each		6341 <input type="checkbox"/>		\$30	
Tray, Height Adj, Clear with Picture Window	Each	66281 <input type="checkbox"/>		67281 <input type="checkbox"/>	\$545	
Mounting Hardware Tray * Req'd For 6628/6728	Each	66282 <input type="checkbox"/>		67282 <input type="checkbox"/>	Included	
Arm Rest, Height, Depth, Angle Adjustable	Pair	8431 <input type="checkbox"/>		8478 <input type="checkbox"/>	\$525	

Upper Body Positioning

Chest Butterfly Style Neoprene	Each	6385 <input type="checkbox"/>		6585 <input type="checkbox"/>	\$215	
Chest Vest Style	Each	6631 <input type="checkbox"/>	6632 <input type="checkbox"/>	6731 <input type="checkbox"/>	\$185	
Chest H-Style, Posture Support, Bodypoint®	Each		214806 <input type="checkbox"/>		\$225	
PivotFit Shoulder Harness, Dynamic, Bodypoint®	Each		214270 <input type="checkbox"/>		\$225	
PivotFit Shoulder Harness, Bodypoint®	Each		214250 <input type="checkbox"/>		\$200	
Trunk Support Lateral Dynamic MIS Microstim	Pair	6383 <input type="checkbox"/>		6583 <input type="checkbox"/>	\$360	
Trunk Support, Rigid Lateral	Pair	6677 <input type="checkbox"/>		6577 <input type="checkbox"/>	\$285	
Trunk Support, Flexible Lateral with Chest Strap	Pair		6596 <input type="checkbox"/>		\$400	

General Accessories

Grab Bar, Padded Black	Each	6621 <input type="checkbox"/>		6721 <input type="checkbox"/>	\$140	
Vent Tray	Each		110004 <input type="checkbox"/>		\$525	
IV Hanger With Mounting Clamp	Each		665011 <input type="checkbox"/>		\$200	
Oxygen Tank Holder, E-Tank	Each		665012 <input type="checkbox"/>		\$200	
Straps, Secure folded seat/frame	Each		667011 <input type="checkbox"/>		\$35	
Straps, Wrist "Runaway" Tether	Each		667111 <input type="checkbox"/>		\$35	

Weather Protection

Canopy, Sun/Rain with Window Black	Each		663761 <input type="checkbox"/>		\$375	
Lower Body Wind Breaker, black for 663761	Each		661852 <input type="checkbox"/>		\$130	
Canopy Sun/Rain with Window Grey	Each		637811 <input type="checkbox"/>		\$415	
Lower Body Wind Breaker, grey for 637811	Each		637812 <input type="checkbox"/>		\$135	
Canopy, Rain Only, Clear	Each		6638 <input type="checkbox"/>		\$250	
Rain Shield, Clear for Canopy 663761/637811	Each		6739 <input type="checkbox"/>		\$150	
Rain cape, Grey	Each	6340 <input type="checkbox"/>		6540 <input type="checkbox"/>	\$200	
Sleeping bag, Summer Weight, Black	Each	6619 <input type="checkbox"/>		6719 <input type="checkbox"/>	\$385	
Sleeping Bag, Winter Weight, Black	Each	6618 <input type="checkbox"/>		6718 <input type="checkbox"/>	\$400	
Umbrella, Angle Adjustable	Each		6594 <input type="checkbox"/>		\$115	