



RAIZER[®]

by Liftup

State of the Art Lifting Chair

Liftup's patented emergency lifting chair is a fast and safe solution when a person has fallen and needs help with transfer to either a standing or sitting position.

Equal Opportunities

Every day thousands of elderly and physically challenged people are falling without being able to get up without assistance. The Raizer from Liftup can with ease be assembled and operated by only one assistant. The fallen person can be raised in a safe and effortless manner within a few minutes.

Easy Operation

The assistant can safely operate the Raizer with minimal physical effort while at the same time observing and caring for the fallen by for instance supporting the head of the fallen person when the transfer is in progress. The Raizer is an optimal solution for personnel in home care as well as for ambulance services and all personnel working with handling and transfer of individuals in general. The Raizer is light and easy to transport. It is carried disassembled to the place where it is to be used. The Raizer is easily assembled and placed around the

fallen individual without strain to neither the fallen person nor the assistant. The Raizer is battery operated and has operation buttons on the side of the chair – and it can be supplied with a remote control as supplementary equipment.

Safety above all

For a fallen individual it is of utmost importance that there is a feeling of security while the transfer to standing or sitting position is taking place. At all times during the transfer with Raizer the individual has contact with the floor as well as physical contact and eye contact with the assistant. For the assistant a transfer can be done without any strain to his or her back. The Raizer can transfer individuals weighing up to 150 kgs. The transfer movement is smooth and protective and as an extra comfort the Raizer is equipped with a seat belt.

BENEFITS

- Mobile chair for transfer from floor to a sitting or standing position.
- Minimal physical effort for the assistant
- Both the fallen person and the assistant can feel comfortable and safe in the process
- Easy to transport in 2 light bags
- Assembly prior to operation: max 3 minutes
- Easy to clean



Technical Specifications

- Lifting capacity: Max. 150 kg
- Lifting time: 20-30 secs.
- Weight:
Seat 9 kg. + legs and back rest 4 kg
– total 13 kg
- Charger (2 pcs.):
12V. 90-264V, < 20W and
12 V car adapter
- Number of lifts on a full charge:
40 with a max. load
Approx. 100 with an average load
- Charging time: Max. 6 hours
- Charging of empty battery:
10–15 min = 1 transfer
- CE-mark
- Patent pending
- Registered design

RAIZER



SCAN
for further info &
Find nearest Distributor

Extras for Raizer:

- Remote control
- Support Belt for Torso
- Hygiene Cover
- Trolley
- Stair Wheels (for trolley)
- Bike bracket (for trolley)

Distributed in Canada by Lighthouse Medical Ltd.
info@lighthousemedicaltd.com or 1.866.609.6827