





Date:	Contact:	PO#	Marked For:
Customer Name:			
Address:			
Confirmation Email To (Buyer):			

Specifications	30400PLEF	30500PLHE	30500SAE
<p>Protekt Patient Lifters</p>  <p>All Protekt Lifts come with Hi-capacity batteries(24V) which last 2x as standard batteries.</p>	<p>Full Body Patient Lift Compact</p>  <p>Qty <input type="text"/></p> <p>The Proactive full body folding patient lift represents a new standard in full body lifts for individuals who travel or have storage issues in their home. The small foot print is ideal for home use. It easily folds up to fit in most automobile trunks for transport or travel.</p>	<p>Full Body Patient Lift</p>  <p>Qty <input type="text"/></p> <p>The Protekt® 500 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.</p>	<p>Sit to Stand Patient Lift</p>  <p>Qty <input type="text"/></p> <p>This 500lb Sit-to-Stand patient lift has a higher weight capacity than most other competitors sit-to-stand lifts. With its 3" dual front casters it ensures safe and easy transfers. Our performance class actuator also includes an emergency quick release.</p>
PRICE	\$3995	\$4995	\$5295
Max Capacity	400lb/180kg	500lb/225kg	500lb/225kg
Base Legs Length	41.6"	44.9"	40"
Base Width - Closed	25.8"	24.9"	26"
Base Width - Open	37.7"	42.9"	38"
Mast Height	47.5"	48.3"	NA
Base Height	4.5"	4.5"	NA
Min. Lifting Height	27.4"	25.8"	33"
Max. Lifting Height	67.2"	69.1"	68"

Patient Lift Slings	LHM Code	WT. Cap	Price	Order
Full Body/Mesh - Med - 41.5" x 54.7"	30110	600lbs	\$425	
Full Body/Mesh - Lrg - 45.5" x 60.5"	30111	600lbs	\$425	
Full Body/Solid - Med - 41.5" x 54.7"	30112	600lbs	\$425	
Full Body/Solid - Lrg - 45.5" x 60.5"	30113	600lbs	\$425	
Full Body/Solid - XLrg - 45.5" x 65.3"	30117	600lbs	\$425	
Full Body/Mesh Commode Open - Med - 41.5" x 54.7"	30114	600lbs	\$425	
Full Body/Mesh Commode Open - Lrg - 45.5" x 60.5"	30115	600lbs	\$425	
Full Body/Mesh Commode Open - XLrg - 45.5" x 65.3"	30116	600lbs	\$425	
Divided Leg Sling/Solid - Small	30100P	450lbs	\$445	
Divided Leg Sling/Solid - Medium	30100	450lbs	\$445	
Divided Leg Sling/Solid - Large	30101	450lbs	\$445	
Divided Leg Sling/Solid - X Large	30102	450lbs	\$445	



2 years - electronics
3 years - frame