

# PROTEKT® 500 LIFT

## 33500 | POWER PATIENT LIFT

The high-performance Propekt® 500 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.

- **500 lb. weight capacity.**
- **6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Capable of lifting from the floor.**
- **User friendly easy to operate base opening control adjusts and locks for added security.**
- **New upgraded electronic controls.**
- **Manual emergency lowering provides safe response to power loss.**
- **Dual front locking casters.**
- **Clip free sling hookups improve safety and reduce downtime.**



(Optional Scale Item # 30300)

### SPECIFICATION FOR ITEM # 33500

Weight Capacity	500 lbs.
Base Legs Length	44.9"
Base Closed External Width	24.9"
Base Open External Width	42.9"
Mast Height	48.3"
Base Height	4.5"
Min. Lifting Height	25.8"
Max. Lifting Height	69.1"
Net Weight	94 lbs.
Shipping Weight	106 lbs.

**WARRANTY:**  
Electronics 2 years  
Frame 3 years



- **Modular design allows easy parts replacement.**
- **Advanced electronics with auxiliary up/down lift controls.**